



**Cementation recommendations for IPS Empress**

Resin cements		Empress Esthetic	e.max full contour and e.max crowns
Product Name	Manufacturer	Veneers	
RelyX ARC	3M Espe	No	Yes
RelyX Veneer	3M Espe	Yes	Yes
C&B Cement	Bisco	Yes	Yes
Choice	Bisco	Yes	Yes
Illusion	Bisco	Yes	Yes
Prevue	Cosmodent	Yes	Yes
Insure	Cosmodent	Yes	Yes
Ultra-Bond Quick	Denmat	Yes	Yes
Calibra	Dentsply	Yes	Yes
Appeal	Ivoclar Vivadent	Yes	Yes
Multilink	Ivoclar Vivadent	No	Yes
Variolink II	Ivoclar Vivadent	Yes	Yes
Nexus	Kerr	Yes	Yes
Panavia 21	Kuraray	No	Yes
Panavia F	Kuraray	No	Yes
C&M Metabond	Parkell	Yes	Yes
Cement-it	Pentron	Yes	Yes
Lute-it	Pentron	Yes	Yes

**Conventional cements for use with e.max.**

VivaGlass CEM	Ivoclar Vivadent	No	Yes
Fuji 1	GC	No	Yes
Ketsic CEM	3M Espe	No	Yes
RelyX Unicem	3M Espe	No	Yes

**Resin reinforced or Resin Modified Hybrid Glass Ionomer Cements**

**Should not be used with Empress products due to expansion**

RelyX Luting Cement	3M Espe	No	No
Fuji CEM	GC	No	No
Fuji Plus	GC	No	No

\* Are Bondable or Cementable